

Safe and easy stent removal anywhere, any time

Isiris α is a patented device specifically designed for cystoscopic ureteral stent removal.

Isiris α is a single use flexible cystoscope, with integrated grasper which connects to a portable monitor to produce digital quality visualization and image recording. *Isiris α* provides clinicians with the choice and opportunity to remove JJ stents in almost any place at any time.

With *Isiris α*, stents can be removed by a solo operator, releasing time in other busy endoscopic departments. *Isiris α* is sterile and ready to use and has no re-processing or maintenance costs.



References

- Robert Koch Institut recommendations for flexible endoscope reprocessing and accessories : http://www.rki.de/DE/Content/Infekt/Krankenhaushygiene/Kommission/Downloads/Medprod_Rili_2012.pdf?__blob=publicationFile. Web accessed on April 26, 2016
- Malavaud S, Boiteux JP, Coloby P, Bugel H, Verine JL, Conquy S, Doublet JD, Bruyère F. [Flexible cystoscopes: disinfection and microbiological surveillance practices among French urologists]. [Article in French] Prog Urol. 2012 Oct;22(12):731-5.
- A multicentric evaluation of a new single use solution for ureteral stent removal, Traxer O., Giusti G., Patterson J., Palmero J.L., Straub M., De La Rosette J., Abstract session EAU congress Munich 2016.
- N° DGOS/PF2/2014/338 du 8 décembre 2014 : French national audit on flexible, non-steam sterilisable endoscopes with a working channel. [<http://www.cclin-arlin.fr/nosobase/Reglementation/2014/Instruction/08122014.pdf> – Web accessed on April 26, 2016]
- ECRI Institute. Endoscope Reprocessing Systems. Tech IQ. 2014.




ISIRIS™ description	Brand name	REF
Single use flexible cystoscope with an integrated grasper (pack of 5)	Isiris™α	ALFA01
Monitor	Isiris™ Monitor	MN0001
Spare part description		
Bracket (Hook+key)	NA	BR0001

PORGÈS

 **Coloplast**

The Coloplast logo and the Porgès logo are registered trademarks of Coloplast A/S. © [2016-06]
All rights reserved. Coloplast A/S, 3050 Humlebaek, Denmark.

Global marketing: DOC0053EN



Isiris™

The single use stent
removal solution




 Coloplast

ISIRIS
PERFORMANCE
VERSATILITY
OPTIMISATION
SAFETY

ISIRIS α is an innovative digital solution to remove ureteral stents. It offers multiple clinical choices to surgeons, reduces steps and simplifies the procedure. ISIRIS α is a patented single use device which is easy to use.

SINGLE USE CYSTOSCOPE

STENT REMOVAL



IsirisTM THE



Revolutionary

PERFORMANCE

- Integrated solution
- State of the art digital technology
- Procedure specific




Freedom to operate

VERSATILITY

- Anywhere anytime
- Solo operator
- Easy to transport and set up

al stents anywhere and anytime. *ISIRIS α* provides greater
the cumbersome environment for stent removal.
use and transport.



MOVAL SOLUTION

MONITOR



Plug & play



DIGITAL IMAGE

PROCEDURE RECORDING

EASY USER INTERFACE

PORTABILITY



THE 4 PILLARS



Resource saving

- Staff time
- Reprocessing and maintenance
- Operating capacity

OPTIMISATION



Protection

- Patient comfort
- No cross contamination
- Equipment preservation from grasper damage

SAFETY